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Cobo A., Outcomes of vitrified early cleavage-stage and blastocyst-stage embryos in a cryopreservation program: evaluation of 3,150 warming cycles. *Fertility & Sterility*, 2012.

Cobo A., Outcome of cryotransfer of embryos developed from vitrified oocytes: double vitrification has no impact on delivery rates. *Fertility & Sterility*, 2013.

Yang H., Comparison of differences in development potentials between frozen-thawed D5 and D6 blastocysts and their relationship with pregnancy outcomes. *Journal Assisted Reproduction Genetics*, 2016.

Coello A., Analysis of the morphological dynamics of blastocysts after vitrification/warming: defining new predictive variables of implantation. *Fertility & Sterility*, 2017.

Liu H., Elevated incidence of monozygotic twinning is associated with extended embryo culture, but not with zona manipulation or freeze-thaw procedure. *Fertility & Sterility*, 2018.

Gu Fang., Perinatal outcomes of singletons following vitrification versus slow-freezing of embryos: a multicenter cohort study using propensity score analysis. *Human Reproduction*, 2019.

La Marca A., A novel transnational fresh oocyte donation (TOD) program based on transport of frozen sperm and embryos. *Human Reproduction*, 2019.

Liu H., Effect of endometrial thickness on ectopic pregnancy in frozen embryo transfer cycles: an analysis including 17,244 pregnancy cycles. *Fertility & Sterility*, 2019.

Zhang J., The impact of embryo quality on singleton birthweight in vitrified-thawed single blastocyst transfer cycles. *Human reproduction*, 2020.

Cobo A., Viral screening of spent culture media and liquid nitrogen samples of oocytes and embryos from hepatitis B, hepatitis C, and human immunodeficiency virus chronically infected women undergoing in vitro fertilization cycles. *Fertility & Sterility*, 2012.

Li W., Influence of storage time on vitrified human cleavage-stage embryos frozen in open system. *Gynecological endocrinology*, 2017.

Sekhon L., Blastocyst vitrification, cryostorage and warming does not affect live birth rate, infant birth weight or timing of delivery. *Reproductive Biomedicine Online*, 2018.

Ueno S., Cryostorage duration does not affect pregnancy and neonatal outcomes: a retrospective single-centre cohort study of vitrified-warmed blastocysts. *Reproductive Biomedicine Online*, 2018.

Costa-Borges N., Cryopreservation of oocytes and embryos in times of COVID-19: Can the cure be worse than the disease? Risk assessment in the IVF laboratory. *Fertility and Sterility*, 2020.

Coello A., Prediction of embryo survival and live birth rates after cryotransfers of vitrified blastocysts, *RBMO*, 2021.

Coticchio G., Perturbations of morphogenesis at the compaction stage affect blastocyst implantation and live birth rates, *Human Reproduction*, 2021.

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Cobo A., Viral screening of spent culture media and liquid nitrogen samples of oocytes and embryos from hepatitis B, hepatitis C, and human immunodeficiency virus chronically infected women undergoing in vitro fertilization cycles. *Fertility & Sterility*, 2012.

Cobo A., Oocyte vitrification as an efficient option for elective fertility preservation. *Fertility & Sterility*, 2016.

Cobo A., Effect of oocyte vitrification on embryo quality: time-lapse analysis and morphokinetic evaluation. *Fertility & Sterility*, 2017.

Diaz-García C., Oocyte vitrification versus ovarian cortex transplantation in fertility preservation for adult women undergoing gonadotoxic treatments: a prospective cohort study. *Fertility & Sterility*, 2018.

Grynberg M., BRCA1/2 gene mutations do not affect the capacity of oocytes from breast cancer candidates for fertility preservation to mature in vitro. *Human Reproduction*, 2018.

Cobo A., Elective and Onco-fertility preservation: factors related to IVF outcomes. *Human Reproduction*, 2018.

Creux H., Thirteen years' experience in fertility preservation for cancer patients after in vitro fertilization and in vitro maturation treatments. *Journal Assisted Reproduction Genetics*, 2018.

Coello A., Effect of oocyte morphology on post-warming survival and embryo development in vitrified autologous oocytes. *Reproductive Biomedicine Online*, 2019.

Cobo A., Use of cryo-banked oocytes in an ovum donation program: a prospective, randomized, controlled, clinical trial. *Human Reproduction*, 2010.

Solé M., How does vitrification affect oocyte viability in oocyte donation cycles? A prospective study to compare outcomes achieved with fresh versus vitrified sibling oocytes. *Human Reproduction*, 2013.

Bárcena P., Should we worry about the clock? Relationship between time to ICSI and reproductive outcomes in cycles with fresh and vitrified oocytes. *Human Reproduction*, 2016.

Domingues TS., Egg donation of vitrified oocytes bank produces similar pregnancy rates by blastocyst transfer when compared to fresh cycle. *Journal Assisted Reproduction Genetics*, 2017.

Parmegiani L., Transnational oocyte donation program: fresh versus vitrified oocytes. *Human Reproduction*, 2019 (Reply Letter).

Cornet-Bartolomé D., Efficiency and efficacy of vitrification in 35.654 sibling oocytes from donation cycles, *Human Reproduction*, 2020.

Costa-Borges N., Cryopreservation of oocytes and embryos in times of COVID-19: Can the cure be worse than the disease? Risk assessment in the IVF laboratory. *Fertility and Sterility*, 2020.

Rienzi L., Definition of a clinical strategy to enhance the efficacy, efficiency and safety of egg donation cycles with imported vitrified oocytes. *Human Reproduction*, 2020.

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Almeida Ferreira Braga DP., Freeze-all, oocyte vitrification, or fresh embryo transfer? Lessons from an egg-sharing donation program. *Fertility & Sterility*, 2016.

Berkkanoglu M., Optimal embryo transfer strategy in poor response may include freeze-all. *Journal Assisted Reproduction Genetics*, 2017.

Xue Y., Freeze-all embryo strategy in poor ovarian responders undergoing ovarian stimulation for in vitro fertilization. *Gynecological Endocrinology*, 2018.

Cardenas Arma D.F., Frozen-thawed blastocyst transfer in natural cycle increase implantation rates compared artificial cycle. *Gynecological Endocrinology*, 2019.

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Ubaldi F.M., Reduction of multiple pregnancies in the advanced maternal age population after implementation of an elective single embryo transfer policy coupled with enhanced embryo selection: pre- and post-intervention study. *Human Reproduction*, 2015.

Rodríguez-Purata J., Reproductive outcome is optimized by genomic embryo screening, vitrification, and subsequent transfer into a prepared synchronous endometrium. *Journal Assisted Reproduction Genetics*, 2016.

Chamayou S., The Accumulation of Vitrified Oocytes Is a Strategy to Increase the Number of Euploid Available Blastocysts for Transfer After Preimplantation Genetic Testing. *Journal Assisted Reproduction Genetics*, 2017.

Cimadomo D., Associations of blastocyst features, trophectoderm biopsy and other laboratory practice with post-warming behavior and implantation. *Human Reproduction*, 2018.

Coll L., Transition from blastomere to trophectoderm biopsy: comparing two preimplantation genetic testing for aneuploidies strategies. *Zygote*, 2018.

Hernandez-Nieto C., What is the reproductive potential of day 7 euploid embryos?. *Human Reproduction*, 2019.

Magli MC., Deoxyribonucleic acid detection in blastocoelic fluid: a new predictor of embryo ploidy and viable pregnancy. *Fertility & Sterility*, 2019.

Wilding M., Thaw, biopsy and refreeze strategy for PGT-A on previously cryopreserved embryos, *Facts Views Vis Obgyn*, 2019.

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